

5.0 8/A herda 10/17/00

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Applicant: Katia Georgopoulos et al.

Art Unit: 1632

Serial No.:

09/259,389

Examiner: Jospeh Woitach

Filed Title February 26, 1999 THE HELIOS GENE

Commissioner for Patents Washington, D.C. 20231

RESPONSE

This is in response to the action dated March 29, 2000. An appropriate extension of time is being filed concurrently herewith. Please amend the application as follows:

In the Specification:

At page 56, line 30, after "Helios polypeptide." please add -As indicated in the cited text from *Current Protocols in Molecular Biology*, a hybridization solution ("hybridization solution I") can, for example, contain (for a total volume of 1 liter): 480 ml formamide; 240 ml 20x SSC; 10 ml 2 M Tris.Cl, pH 7.6; 10 ml 100x Denhardts solution; 50 ml deionized water; 200 ml 50% dextran sulfate; and 10 ml 10% SDS. A low stringency wash buffer that can be used after hybridization in hybridization solution I can contain: 2x SSC; and 0.1% SDS. A high stringency wash buffer that could be used after hybridization in hybridization solution I can contain: 0.2x SSC; and 0.1% sodium dodecyl sulfate (SDS).

An alternative hybridization solution ("hybridization solution II") can contain: 1% crystalline bovine serum albumin (BSA) (fraction V); 1 mM EDTA; 0.5 M NaHP0₄, pH 7.2; and 7% SDS. A low stringency wash buffer that can be used after hybridization in hybridization solution II can, for example, contain: 0.5% BSA (fraction V); 1 mM Na₂EDTA; 40 mM NaHPO₄, pH 7.2; and 5% SDS. A high stringency wash buffer that can be used after

CERTIFICATE OF MAILING BY FIRST CLASS MAIL

I hereby certify under 37 CFR §1.8(a) that this correspondence is being deposited with the United States Postal Service as first class mail with sufficient postage on the date indicated below and is addressed to the Commissioner for Patents, Washington, D.C. 20231

1 / 0 d / 0 map 011

Signature

Typed or Printed Name of Person Signing Certificate

(h)